Product Information

CONFORM® Med. Air Cushions are made of natural rubber. Natural Rubber is a renewable resource which requires only low energy expenditure processes for its treatment. CONFORM® Products are fully biodegradable and require no special disposal measure.

Only the purest rubber (24-015 NR green) is used for the production of CONFORM® products. This particular quality is free of allergens and fulfils all current demands for a medical product.

CONFORM® cushions are available in 3 different heights (5, 7.5 and 10 cm) and a wide range of sizes. Our cushions offer the very highest seating comfort, limiting the danger of pressure sores, and can be used therapeutically with existing skin problems. The high elasticity results in pressure relief - allowing the blood flow in the capillaries. The cushions are used from stage of prophylaxis up to pressure ulcer (Grade IV).

There are no usage restrictions either in terms of the material or the cushion form. The choice of 2- or 4-valve cushions allows the seating comfort to be optimized when necessary.

Inaccuracy in the application and use of the product and consequent risks remain those of the "user". No contraindications have been reported till date.

Maximum inflation pressure (cushion only) 20-30mb.

Maximum users weight - no limitation. For patients (weight above 190 kg) we recommend at least 7.5 or 10 cm height cushion.

The cushion is equipped with a flame retardant cotton/lycra-stretch cover with non-slip base. The cushion corners are designed with integrated metal eyelets to allow for fixing strips etc.

Cushion weights vary depending on size i.e.:

Code	Measurement/Height	Weight
SKN4040 (1valve)	38 x 38 x 5 cm	850 g
SKN4045 (1valve)	38 x 43 x 5 cm	1020 g
SKD4540 (2valve)	43 x 38 x 7.5 cm	1040 g
SKD5045 (2valve)	47 x 43 x 7.5 cm	1600 g
SKH3540 (4valve)	33 x 38 x 10 cm	1090 g
SKH5055 (4valve)	47 x 52 x 10 cm	2230 g

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CONFORM® Air Cushions are Medical Products class I (DE/CA/23/051603NA001.01)

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CE-Declaration:

CONFORM® Air Cushion Systems are in accordance with all requirements for medical products as set out in 93/42/EWG, made in the EC.

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DIN EN ISO 9001:2008
DIN EN ISO 13485:2010

**PRODUCT INFORMATION
OPERATING INSTRUCTIONS**

Operating Instructions

Inflate CONFORM® Air Cushion with the enclosed pump until it begins to arch. Attach the pump tubing to the cushion valve. Open the valve by turning it anticlockwise and inflate the cushion. Close the valve by turning it clockwise.

Seat the patient and deflate the cushion so that the patient comfortably sinks, but without allowing the patients bottom to touch the base of the cushion.

The distance between the deepest point of a sitting person (usually the ischial tuberosities) and the base of cushion should be approximately ¾ of an inch or 2 cm.

If it is possible to insert your hand between the patients bottom and the base of the cushion, it has been set correctly.

Cleaning and maintenance

The use of only CONFORM® covers with the CONFORM® cushion is recommended. Natural rubber should not be exposed to sunlight (UV-light), which could result in premature ageing of the material. Please do not through the case away. When your CONFORM® cushion is not in use, it's box is required for storage.

Both cover and cushion can be wash from 40 – 60° Celsius. We recommend machine washing from 30°C - 40°C (without spin cycle) for the cushion. Please ensure that the valves of the cushion are closed during this process.

CONFORM® Air Cushions are salt-water resistant.

The requirement for storage, cleaning and maintenance are in accordance with ISO 2230.

Warranty

CONFORM® Air Cushion Systems have a manufacturing warranty of 3 years.

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