



Pro-care Auto

High Specification Pressure Relieving Mattress

APEX
Care for a Healthy Life

Pro-care Auto:

High Specification Pressure Relieving Dynamic Mattress

The Pro-care Auto is a high specification pressure relieving dynamic mattress designed to offer exceptional levels of support in any healthcare environment.

The Pro-care Auto dynamic mattress system automatically calibrates in accordance to the patients weight and position in the bed, eliminating the risk of injury to the patient by the caregiver entering the patients weight incorrectly.

The Pro-care Auto has a range of unique features for both patients and caregivers, and has been designed to comply with the most stringent infection control guidelines. It is suitable for all types of patient from pressure ulcer prevention right the way through to managing the most severe pressure ulcers.

The Pro-care Auto is one of the only mattresses on the market that has a designated heel zone allowing the caregiver to float the heel, meaning they can comply with EPUAP guidelines on pressure ulcer management in the heel area.



CPR: The CPR dial is located at the head end of the mattress and is designed to be highly visible in an emergency. Its unique dial design reduces the risk of accidental activation.



Transport Mode: Transport Mode is designed to be simple and easy to use. Its simple click in design also means it is robust and easy to keep clean.



Antimicrobial Filter: The antimicrobial filter is designed to ensure the air going into the system is as clean as possible, helping to protect the pump from break down and to ensure air stored in the mattress is as clean as possible.



Air Cells: The air cells modular design means maintenance, replacement and decontamination are easy to carry out. Cell on Cell Technology also means that the patient is safely supported throughout the length of the mattress, and the unique rectangular cell design in the torso section means the patient benefits from higher levels of pressure relief and weight redistribution.

Product features and benefits

Designated Heel Zone: Allows the caregiver to isolate individual air cells in the heel area for 100% pressure offloading, allowing the caregiver to float the heel and comply with EPUAP guidelines relating to pressure ulcer management.

Auto Adjust: Automatically adjusts the internal cell pressures depending on the individual patients weight and position on the mattress, eliminating the need for the caregiver to manually set the pump to the patients weight.

Alternating Therapeutic Mode: Continuously and sequentially inflates and deflates the air cells to allow for pressure redistribution over a ten minute cycle time.

Static Mode: Addresses areas of high contact pressure by redistributing the body mass over a greater surface area through a process of immersion. The mattress will not stay in the static mode but will revert back to the therapeutic mode after 20 minutes.

Seat Inflate Mode: When a patient is in the upright position, the seat inflate mode increases the pressure throughout the mattress in order to better support the sacral area.

Low Pressure Alarm: Alerts the caregiver if there is a loss of air pressure in the mattress.

Max Firm Mode: Hyper inflates the cells providing a firm surface for moving, handling and nursing procedures. For the safety of the patient the mattress will only stay in Max Firm Mode for 20 minutes and will then revert back to Alternating Therapeutic Mode.

Power Failure Alarm: Notifies the caregiver when there is a loss of power if the mattress is unplugged.

Service Fault Indicator: Alerts the caregiver when the mattress needs to be repaired.

Lock Out: Automatically locks the control panel to help prevent the changing of settings inappropriately.

Rectangular Sacral Cell on Cell Design: Unique cell design offers the patient more support and comfort in the Sacral area throughout the therapeutic pressure relieving cycle.


Mattress Cover: The cover is flame resistant to BS 7175 Crib 7 and has a host of features to aid with infection prevention including a 360 degree zip, white underside, vapour permeability, antimicrobial properties and is water resistant. The cover is four way stretch to offer greater pressure relief and comfort levels to the patient, it is also possible to launder the cover in a washing machine.

Mattress Specification	
Mattress Dimensions	: 200x85x20.3 cm or 200x90x20.3 cm
Mattress Weight	: 7.1 KG
CPR	: Robust CPR Valve
Cell Structure	: Cell on Cell Technology
Designated Heel Zone	: Allows for complete offloading/floating of heels
Max Firm	: To assist with patient care
Seat Inflation Mode	: Helps support the patient in the upright position
Cable Management	: To allow for safe storage of the power cable
Modular Cell Construction	: For ease of maintenance
Anchoring Straps	: For safe attachment to bed frame
Static Head Zone	
Anti-Sink Rectangular Sacral Support Cells	
Static Mode	
Non-Slip Base	


Cover Specification	
Flame Resistant Crib 7	
360 Degree Zip	
White Underside of Cover	
Four Way Stretch	
Vapour Permeable	
Water Resistant	
Machine Washable	
Antimicrobial	

Pump Specification	
Maximum Patient Weight	: 250 KG
Pump Size	: 29.1x20x11.7 cm
Pump Weight	: 2.2 KG
Cycle Time	: 10 minutes
Comfort Control	: Manual comfort level adjustment
Voltage	: AC220-240V, 50Hz; AC : 120V/60Hz
ABS Cover	: ABS Flame Retardant Pump Cover
Transport Mode	
Hanging Brackets	
Audible and Visual System Fault Alarms	
Antimicrobial Air Filter	
Extra Long Power Cable	
Mute Button	
Easy to Identify Power Light	
Pump Lockout	
Fault Indicator	

Alternating Support Surface



Pressure Ulcer Classification Guide



Authorised Distributer

O'Flynn Medical

The Old Creamery, Tanyard,
Millstreet Town,
Co Cork. P51 D328.
T: 029 21799
Lo Call: 1890 440 440
E: info@oflynnmedical.com



www.oflynnmedical.com