

Your residents are unique.
Are your fall management solutions?



O'FLYNN MEDICAL Ltd

1890 440 440

www.oflynnmedical.com



Each one of your residents is unique. They each have different capabilities and require different levels of care and attention. We designed the CareSense line of fall management systems to provide a high level of flexibility with numerous monitor and sensor options so you can create a customized solution that's best for each resident.

Meet Frank

Frank frequently gets out of bed at night. The staff needs a solution that notifies them of potential exits without disturbing nearby residents.



All CareSense monitors feature a jack for connection to a nurse call system, so a potential fall can be signalled at the nurse's station without a loud alarm. The Advanced Bed and Chair Monitors (CSM-BC200 and CSM-BC250) also feature a Volume button to quickly change the alarm volume, and a Reset button to quickly cancel an alarm.

Meet Joe

Joe has tripped over cords in the past. The staff at his facility needs a solution to keep this from happening again.



When used with wireless sensor pads the Wireless Advanced LCD Bed and Chair Monitor (CSM-BC400) eliminates cord tripping hazards, as well as cord damage. With its easy-to-read LCD screen, dedicated function buttons, and wireless capabilities, it's the ideal choice when you need a fall management monitor that's flexible enough to use with any of your residents in any situation.

Bed and Chair Sensor Pads

To meet the needs of different residents and situations, CareSense pads are available in different sizes and styles.

- ▶ Corded pads feature a breakaway cable
- ▶ Timed pads signal the monitor when their usable life is nearing its end
- ▶ Wireless pads eliminate tripping hazards and cord damage



20" x 30" Bed Pad



10" x 15" Chair Pad

Bed and Chair Pad Styles

Corded Pads

CareSense corded pads are a cost-effective choice when you need a simple, familiar sensor for your monitors.

Wireless Pads

Wireless pads (for use with the Wireless Advanced LCD Monitor) eliminate cord tripping hazards and cord damage — the primary cause of false and failed alarms.

	20" x 30" Bed Pads	10" x 15" Chair Pads
	1 Year	1 Year
Non-Timed Corded Pad	●	●
Timed Wireless Pad	●	●

Additional Sensors



Seat Belts

CareSense Seat Belts are available in buckle and Velcro® versions. In addition to setting off an alarm when unfastened, they also help keep residents from sliding down in their chairs.



Toilet Sensor Pad



The CareSense Toilet Sensor Pad mounts discreetly underneath the toilet seat so that residents can be monitored by staff without a loss of privacy.



Floor Mats and Pads

CareSense Floor Mats can be placed near beds, chairs, and doorways to alert staff if a resident attempts to stand or walk. Floor pads provide extra protection with 2" padding

Features and Compatibility

	Advanced 200 (CSM-BC200)	Wireless Advanced (CSM-BC400)
Monitor Features		
Two-Year Warranty	●	●
Wall Mounting Holes and Hardware	●	●
Cord Strain Relief	●	●
Clip and Velcro Hanger	●	●
Nurse Call System Interface	●	●
Call Cord Receptacle	●	●
Battery Powered	●	●
AC Powered		
Optional AC Power Adaptor	●	●
Battery Back Up		
On/Off Switch		
Reset Button	●	●
Volume Button	●	●
Adjustable Volume	●	●
Selectable Tone Settings		●
Low Timed Pad Notification	●	●
Low Battery Notification	●	●
Adjustable Time Delay	●	●
Tamper Resistant Option (CareSense Key)	●	●
30 Second Pause Option	●	●
In Use LED Indicator	●	●
Alarm/Low Battery LED	●	●
Pad/Signal Lost LED	●	●
Side-Mounted LEDs		●
LCD Display		●

Sensor Compatibility

Standard Corded Pads	●	
Timed Corded Pads	●	●
Wireless Pads		●
Corded Floor Mat	●	●
Wireless Floor Mat		●
Floor Cushion Sensor	●	●
Buckle Seat Belt	●	●
Velcro Seat Belt	●	●
Toilet Seat Sensor	●	●