

Eclips XXL O'FLYNN MEDICAL

XXL WHEELCHAIR

Aluminum foldable wheelchair with double cross. The Eclips XXL provides a natural comfort with its anatomical seat and backrest. This wheelchair, equipped with reinforced rear wheels, can support a maximum weight of 200 kg.

Easy transfers thanks to the foldable aluminum footplate in one piece.



C27
blue

200KG



B78 - Anti-tipping




B74 - Brakes



ECLIPS XXL

Wide range of accessories and options.



| | | | | | | | | | | |
|------|-----|-----|-----|-----|-----|------|-----|------|------|-----|
| 850 | 375 | 600 | 200 | 530 | 620 | 1070 | 530 | 1140 | 38,5 | 200 |
| 900 | 375 | 650 | 200 | 530 | 620 | 1070 | 540 | 1140 | 39 | 200 |
| 950 | 375 | 700 | 200 | 530 | 620 | 1070 | 550 | 1140 | 40 | 200 |
| 1000 | 375 | 750 | 200 | 530 | 620 | 1070 | 560 | 1140 | 40,5 | 200 |



Sizes in mm, weight in kg. Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible.