



Aquatec Ocean VIP

**German Innovation Only in Cars?
Think again!**

The AQUATEC OCEAN VIP (Tilt in Space) is a high quality, German-designed stainless steel mobile shower commode with 0-35 degree tilt for the more passive/disabled user. High quality castors and a directional lock ensure a smooth ride for all – even on carpet!



Configurations

- Comes with soft seat, pan, and pan holder
- Additional Backrest supplied Free of Charge
- Warranty - 5 years on frame and up to two years on all other components*

*Soft seat and backrest cover are warranted for six months, all other components (excluding frame) are warranted for a period of 2 years



Seat tilt can be individually adjusted (tilt angle 0° - 35°). Exceptionally stable, even if tilted to maximum position.



Ergonomically-designed, height adjustable seat is removable and easy to clean.

Features	Benefits
High quality stainless steel frame and components	Offering strength and stability when it matters most
High quality, rust-proof castors, with directional lock	Providing durability and ease of use
Unprecedented adjustability – height adjustable seat, armrests and footrests	Making the chair comfortable and safe for both patient and carer
Washable, tension-adjustable back-rest with quick-release buckles	Providing optimum comfort and positioning for the user

INVACARE® **Aquatec Ocean VIP Accessories**



Anti-Tippers
Item no: 1470833



Chest and Pelvis belt
Chest 1470081
Pelvis 1470082



Splash Guard - For soft seat.
Item no: 16511



Safety Bar - Length 540mm.
Item no: 10167-10



Armrest Lock - locks the armrest in
lower position.
Item no: 1470077



Insert for soft seat.
Item no: 16343

INVACARE BATHROOM AIDS



INVACARE®

Aquatec Ocean

**German Innovation Only in Cars?
Think again!**

With its unprecedented adjustability, and the smoothest castors around (even on carpet!), the AQUATEC OCEAN is a high quality, German-designed stainless steel mobile shower commode with unprecedented adjustability.



OCEAN transit model



OCEAN self propelled model



design award
winner

2004

Configurations

- Comes with soft seat, pan, and pan holder
- Additional Backrest supplied Free of Charge
- Warranty - 5 years on frame and up to two years on all other components*

*Soft seat and backrest cover are warranted for six months, all other components (excluding frame) are warranted for a period of 2 years



Four braked castors for optimum safety.



Ergonomically-designed, height adjustable seat is removable and easy to clean.

Features	Benefits
High quality stainless steel frames and components	Offering strength and stability when it matters most
High quality, rust-proof castors	Providing durability and ease of use
Unprecedented adjustability – height adjustable seat, armrests and footrests	Making the chairs comfortable and safe for both patient and carer
Washable, tension-adjustable back-rest with quick-release buckles	Providing optimum comfort and positioning for the user

INVACARE®

Aquatec Ocean Accessories



Armrest Width Extension (Upgrade Kit OCEAN XL) Increases load capacity to 150kg.
Item no: 16515



Amputation Support (Above Knee) suits left/right side. Adjustable height & angle.
Item no: 16616



Amputation Support - (Below Knee) suits left/right side. Adjustable height & angle.
Item no: 16615



Chest and Pelvis belt
Chest 1470081
Pelvis 1470082



Splash Guard - For soft seat.
Item no: 16511



Safety Bar - Length 540mm. NB: Fits standard OCEAN models only, not XL.
Item no: 10167-10



Armrest Lock locks the armrest in lower position.
Item no: 16597



Insert for soft seat.
Item no: 16343



24" Self-Propelling Wheels (Pair) to retrofit OCEAN Transit.
Item no: 16379

INVACARE®

Aquatec Ocean Range Specifications

Dimensions	AQUATEC OCEAN TRANSIT (7.10.910)	AQUATEC OCEAN SELF-PROPELLED (7.12.910)	AQUATEC OCEAN VIP (1470793)
Seat height	477-600mm	477-600mm	500-600mm
Seat width	480mm	480mm	500mm
Seat depth	450mm	450mm	452mm
Overall width	560mm (640mm XL ext)	675mm	576mm
Overall length	934mm	1060mm	1058mm
Length without footrests	730mm	856mm	848mm
Footrest depth	154mm	154mm	160mm
Overall height	964-1094mm	964-1094mm	1085-1185mm (without head support)
Height to armrests	680-810mm	680-810mm	719-845mm
Height - footrests to seat	280-530mm	280-530mm	360-485mm
Height - ground to lower seat edge with holder for sanitary pan	390-520mm	390-520mm	422-522mm
Height - ground to lower seat edge without holder for sanitary pan	410-540mm	410-540mm	445-545mm
Width between armrests	459mm (539mm XL ext)	459mm (539mm XL ext)	457mm
Width between 610mm (24") wheels	N/A	530mm	N/A
Height adjustment steps	25mm	25mm	50mm
Castor diameter	127mm (5")	127mm (5")	127mm (5")
Tilt angle	N/A	N/A	0° - 35°
Overall weight	14.5kg	19.5kg	21.5kg
Max weight capacity	130kg (150kg XL ext)	130kg (150kg XL ext)	150kg



Yes, you can.™