



NEW

## Ameda® Finesse™ Breast Pump

Hospital technology  
at home



# Ameda Finesse Breast Pump

Hospital technology at home

Ameda has set the standard for breast pump quality and performance for over 75 years; we are the original inventor of the technology still used in hospital breast pumps today. Ameda's new Finesse personal pump features our hospital-trusted pump technology and is designed for moms at home, at work or on-the-go. Ameda Finesse is ideal for those who want a fast, convenient and easy breast pumping experience.

## Finesse Features:

- Lightweight, portable, and can operate on batteries for easy transport
- Easy to assemble and clean
- No need to clean or replace tubing
- Our quietest personal electric pump for a more discreet experience
- Long, reliable performance backed by a 2-year warranty

## Achieve more comfort

Independent speed and suction controls allow you to choose settings most comfortable for an increased milk supply.



## Keep milk safe

Some leading breast pumps don't prevent moisture and milk from getting into the tubing. Ameda's true closed-system pump kit has the only FDA-cleared barrier that helps protect breast milk and baby from bacteria and viruses while pumping. Without an effective barrier, mold can grow in the tubing and pump.



### COMPETING BRAND



Moisture in Tube = Mold



No Moisture in Tube = No Mold

AMEDA BRAND



## Establish more milk

As used in Ameda's hospital pumps, ComfortFlow™ Technology is designed to mimic baby's natural suckling. This technology provides a comfortable and smooth feel for mom while pumping which is clinically proven to establish a better milk flow.

### Need help choosing a breast pump?

Ameda's ParentCare team is here if you have any questions about Ameda products. Additionally, we have Certified Lactation Counselors (CLCs), an International Board-Certified Lactation Consultant (IBCLC) and Registered Nurse (RN) able to provide lactation guidance and support. **1.866.99AMEDA (1.866.992.6332) or ParentCare@ameda.com**

Qualify for a **FREE** Ameda Pump through your insurance or purchase on Amazon.

For more info: [Ameda.com/Finesse](https://www.ameda.com/Finesse)

Ameda, Inc.  
485 Half Day Road  
Buffalo Grove, IL 60089 USA

1.866.99AMEDA

The Ameda logo, Ameda, Ameda Platinum, ComfortFlow, CustomControl, CustomFit Flange System, Finesse, HygieniKit, and Proven Airlock Protection are trademarks and/or registered trademarks of Ameda, Inc. in the United States and other countries. 1500092-0318

[Ameda.com](https://www.ameda.com) || 1.866.99AMEDA (1.866.992.6332)

