

ICON and ICON-Low adjustable care beds from Nexus DMS



ICON & ICON-Low
End-lift, 4-section adjustable profiling beds

Height range: 40cm to 80cm (ICON)
20cm to 60cm (ICON-Low)

- Safe Working Load: 210kg (33st)
- Protective plate-style headboards in attractive beech
- Full-length timber or steel siderails available
- Auto-regression ("floating back-rest") available
- 20cm extension available
- Hand control features function lock, tilt/anti-tilt and auto-profile (back-rest and knee-break can profile simultaneously to prevent sliding during profiling)
- Linak Electronics & Control System
- Steel mesh mattress platform enhances robust structure, spreads load and allows mattress to breathe, prolonging its life.
- 3 year warranty in line with Nexus DMS Conditions of Sale and website "Useful Information"
- Conforms to BS EN 60601-2-52 standard

Esteem enclosed head/foot boards
Optional enclosed timber head/foot boards offering an aesthetic finish for a more homely feel.

Esteem head/foot boards are available for both ICON and ICON-Low.



Made to match: everything you'll need

Complete range of optional bed add-ons, pressure-redistribution mattresses & accessories, fully compatible with ICON & ICON-Low:

Comfort & pressure-redistribution:

- **Motif mattresses:** Low, Medium, High and Very-high risk options – purpose-designed to fit ICON beds, offering complete peace of mind that the chosen bed/mattress combination will fit properly and work correctly.
- **Optima dynamic air-float system:** Unique reactive mattress system aiding pressure ulcer therapy and pain management by minimising interface pressure and tissue strain across the user's entire body.



Moving & repositioning:

- **Easy-Ex movement aid:** Multi-positional movement aid allowing easier access to & from bed
- **Fowler-position leg angle adjuster:** 6 additional leg-rest angles for added comfort
- **Lifting pole & grab handle:** Overhead aid helping users turn/reposition themselves



Safety:

- **Safety foam bumpers:** For head/foot boards, standard and extended-height side-rails
- **Side-rail height extenders:** Extra protection against falls, as well as maintaining conformity to BS EN 60601-2-52 when a deeper lying surface (e.g. mattress overlay/dynamic system) is used.
- **Crash-Ex safety mat:** Preventing harm to users prone to climbing/falling from bed
- **Ten-Ex total-enclosure bumper system:** Additional padded protection against entrapment or self-harm.

Product design may change without notice in order to improve quality – thank you for your understanding.