



### **Flexibility and comfort – an ideal solution for Acute and Residential care**

The *Invacare Jasmine* Mobile Lifter is a flexible solution that offers a high level of comfort for all clients. A versatile lifter, the **Jasmine** provides safe patient handling up to 32 stone (200kg) weight capacity. With an expansive, wide-opening base and high lifting mechanism, the **Jasmine** offers easy manoeuvrability and positioning for a wide range of lifts. With the latest control box technology, the **Jasmine** offers a 'Soft Start' and 'Soft Stop' function that provides a more comfortable and dignified transfer.



200kg weight limit  
5 year warranty

Trusted Solutions, Passionate People





**Features and Options**

**Flexibility and optimized maneuverability**

The *Invacare Jasmine* is exceptionally mobile and easy to transport. The wide-opening base allows easy access for a wheelchair or armchair, whilst the high lifting mechanism is essential for transfers to and from elevated surfaces. Flexible and with optimized maneuverability, the *Jasmine* covers the needs of both the carer and client.

**Safety first, with an intelligent control box**

The *Jasmine* is fitted with the control box with diagnostic system: an intelligent service lifetime monitoring system that highlights the battery status and service condition of the actuator to ensure safe patient handling.



**Spreader bar options**

The *Jasmine* can be ordered with 2 or 4-point spreader bars. In widths: 450 mm and 550 mm.



**Scale**

A scale is available. Approved according to medical standard EN 45501, certified and sealed.



**Linak® Jumbo™ Care**

The intelligent monitoring system includes a service reminder, emergency operation, audio battery alarm, and battery indicator.

**Technical data**



*Invacare Jasmine* 134 cm



69 cm



105 cm



54 - 192 cm



12 cm



*Invacare Jasmine* 153 cm



200 kg




43 kg



40 cycles

**Frame colours**



Moderate white



Blue dusk

**Warranty**

5 year warranty:  
Frame, Control box,  
Actuator

2 year warranty:  
Hand control

1 year warranty:  
Battery

Invacare Limited  
Pencoed Technology Park -  
Pencoed - Bridgend CF35 5AQ -  
United Kingdom  
Tel: +44 1 656 776 222  
Fax: +44 1656 776 220  
E-mail: [uk@invacare.com](mailto:uk@invacare.com)  
Sales Order E-mail:  
[ordersuk@invacare.com](mailto:ordersuk@invacare.com)  
[www.invacare.co.uk](http://www.invacare.co.uk)

Invacare Ireland Ltd  
Unit 5 - Seatown Business  
Campus  
Seatown Road - SWORDS -  
County Dublin - Ireland  
Tel. +353 1 8107084  
Fax +353 1 8107085  
Email: [Ireland@invacare.com](mailto:Ireland@invacare.com)  
[www.invacare.ie](http://www.invacare.ie)

© 2016 Invacare International Sàrl  
All rights reserved.

All information quoted is believed to be correct at time of print. Invacare reserves the right to alter product specifications without prior consultation.

Jasmine - UK - 01/2016  
Brochure re-order code: 1539003

The *Invacare Jasmine* is part of the 'Safe Patient Handling' range.



safePatienthandling



Yes, you can.®