



5 year warranty

Versatile and safe - an ideal solution for long-term care facilities

The *Invacare Roze* stand assist offers an alternative to a passive lift when a degree of weight bearing ability remains. Ideal for rehabilitation, degenerative conditions and clients whose condition varies from day to day, the **Roze** can provide dignity and a greater degree of independence. With the latest control box technology, the **Roze** offers a 'Soft Start' and 'Soft Stop' function that provides a more comfortable transfer. Combining a range of adjustable features with a lightweight aluminum construction, the **Roze** is equally at home in Residential and Acute Care environments.

Trusted Solutions, Passionate People



Yes, you can.®



Features and Options

A versatile stand assist

The *Invacare Roze* is a versatile stand assist lifter that supports the transfer of individuals who can only bear a minimal amount of their own weight. Ideal for transfers from a seated to seated or seated to standing position, the small footprint of the **Roze** and wide leg opening offers a flexible and comfortable transfer.

Safety first, with an intelligent control box

The **Roze** is fitted with a control box that has a diagnostic system: an intelligent service lifetime monitoring system that highlights the battery status and service condition of the actuator to ensure safe patient handling.



Height-adjustable leg support

The padded leg support can be adjusted for different needs.



Removable foot plate

The foot plate is horizontal and can be removed for standing/walking a patient during rehabilitation.



Control box with diagnostic system

This intelligent monitoring system includes a service reminder, emergency operation, audio battery alarm, and battery indicator.

Technical data



Invacare Roze 94 - 161 cm



105 cm



66 cm



111 cm



138 cm



Invacare Roze 12 cm



100-200 lifts per full charge




53 kg



Max. 200 kg

Frame colours



Moderate white



Blue dusk

Warranty

- 5 year warranty:
 - Frame
 - Control box
 - Actuator
- 2 year warranty:
 - Hand set
- 1 year warranty:
 - Battery

Invacare Limited
 Pencoed Technology Park -
 Pencoed - Bridgend CF35 5AQ -
 United Kingdom
 Tel: +44 1656 776 222
 Fax: +44 1656 776 220
 E-mail: uk@invacare.com
 Sales Order E-mail:
 ordersuk@invacare.com
 www.invacare.co.uk

Invacare Ireland Ltd
 Unit 5 - Seatown Business Campus
 Seatown Road - SWORDS - County
 Dublin - Ireland
 Tel. +353 1 8107084
 Fax +353 1 8107085
 Email: ireland@invacare.com
 www.invacare.ie

© 2016 Invacare International Sàrl
 All rights reserved.

All information quoted is believed to be correct at time of print. Invacare reserves the right to alter product specifications without prior consultation.

Roze - UK - 01/2016 - 1539002



Authorised Distributer
O'Flynn Medical
 The Old Creamery, Tanyard,
 Millstreet Town,
 Co Cork. P51 D328.
 T: 029 21799
 Lo Call: 1890 440 440
 E: info@oflynnmedical.com



Yes, you can.®