

286 E


ROLLATOR

4-wheel foldable rollator. Equipped with a removable basket, plexi table, seat, crutch holder. Parking brakes and height adjustable handles.

Very simple and fast folding system.



C59
champagne



780-970

644

605

590

430

11

120



Sizes in mm, weight in kg. Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible.

Quadri-Light

ROLLATOR

A lightweight, foldable 4-wheel aluminium rollator with parkingbrakes. Ideal for people with (lesser) walking capabilities but who still have the strength and power to use the brakes when they are about to lose their balance.

A smooth plastic seat in combination with a supporting belt when the user wishes to rest.



C29
grey

**optional
transportbag
available**



805-930

700

610

650

290

450

7,25

130



Sizes in mm, weight in kg. Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible.