

# 708 D

## O'FLYNN MEDICAL

Standard & Light manual wheelchairs



### Standard wheelchair with a fixed back

The armrests are foldable and removable and equipped with soft, extra long and comfortable padding. The fabric of the seat and back, together with aluminum legrests, will provide optimal comfort while using this chair.



AMP kit  
(amputated version)



708 Delight S  
lowered version



T30 kit - push version



Double handrail




Removable and foldable armrests with a comfortable padding

Rear wheels with a Quickrelease system

Sturdy armrests with an easy locking system

Small side wheels on the footplates will protect your walls and furniture while moving around

### Wide range of accessories and options.



550	280	390	210	510	430	400-540	430	1100	18	120
580	280	420	210	510	430	400-540	430	1100	18	120
600	280	440	210	510	430	400-540	430	1100	19	120
620	280	460	210	510	430	400-540	430	1100	19	120
640	280	480	210	510	430	400-540	430	1100	19	120
660	280	500	210	510	430	400-540	430	1100	19	120

Sizes in mm, weight in kg. Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible.

