


Lift Chair

Ontario I & II - Montreal



Technical information Lift Chair



	total width	seat width	height armrests	seat height	total height	total length	seat depth	height backrest	weight	users weight
Ontario I	800	465	150	500	1160	840	450	695	49,6	150
Ontario II	845	480	140	500	1145	960	530	670	53,3	150
Montreal	850	485	150	500	1140	830	510	700	54	150

Deviations in sizes and weight (15 mm and 1,5 kg) are possible



Ontario I & II

The Ontario liftchair offers you optimal comfort in your own house. With a simple push on a button you can relax, watch tv or have a little nap in your Ontario liftchair.

The **2 motors** version allows you independent movement of the backrest and the legrest. You can

adjust the legrest without the backrest and vice versa.



Montreal

Make your choice out of different colors and textures for an optimal match with your existing interior. This liftchair can be placed very close to the wall (100 mm) and still recline completely. The footrest inclines first making it easy to watch television without having to lay back completely.



Option:
armrest covers

Colors Ontario



Colors Montreal



Authorised Distributer
O'Flynn Medical
 The Old Creamery, Tanyard,
 Millstreet Town,
 Co Cork. P51 D328.
 T: 029 21799
 Lo Call: 1890 440 440
 E: info@oflynnmedical.com

09/2013



www.oflynnmedical.com



VERMEIREN
 Opening up new horizons in mobility

Subject to technical alterations

R.E.: S.A. Vermeiren