



Sedens 410

Pressure Relieving Seat Cushion

Pressure Relief for a 24 Hour Care Package

The Sedens 410 is a light weight and portable seat cushion designed specifically to be used in chairs along with APEX dynamic mattress system pumps. It's able to extend pressure relief therapy to sit position that provides patient 24 hour pressure area care, allowing caregivers to make versatile care plan for patient to prevent pressure injuries and increase patient's mobility.

Care for a Healthy Life

APEX

Sedens 410

Pressure Relieving Seat Cushion

Providing 24 Hour Pressure Area Care for Patient

- Alternating cushion works with APEX dynamic mattress system pumps.
- Hybrid design with alternating therapy to relieve pressure from patient's buttock whilst reducing the risk of bottoming-out.

Effective seating solution to provide great pressure relief to patient when sitting in chair.

Crib 7 fire retardant, high moisture vapour transmission stretch cover with welded seams to wick moisture away from the surface while preventing any fluid ingress.



6 Air-cell-like one-piece pad with foam encased in to reduce risk of bottoming out.

Non-slip base prevents the cushion from sliding out of position, giving peace of mind to the user.

Specifications

Cushion	Cushion Dimension	44 x 41 x 9 cm
	Cell Height	3"
	Weight	1.11 kg
	Cover Material	Stretch PU Polyester
	Cell Material	Nylon PU
	Base Material	Non-slip PU Polyester
	Fire Retardant Rating	EN 597-1 & EN 597-2
	Maximum Patient Weight	No weight limit, yet the cushion must be properly sized to the individual

Specifications are subject to change for improvement without notice

Apex Medical Corp.

No. 9, Min Sheng St., Tu-Cheng,
New Taipei City, 23679, Taiwan
T +886 2 2268 3100
F +886 2 2268 6525
marketing@apexmedicalcorp.com